THERAPY TIPS

RELIEVE BODY/MUSCLE TENSION

Your Hot and Cold Massage Kit works to relieve tension through muscular massage on your back, shoulders, and glutes. To perform this massage, lay one of the therapy balls on an even, non-slippery surface (e.g. carpet, exercise mat), and roll it beneath the tight or sore body part. To better relieve sore upper back and shoulders, place the ball between the wall and the affected body part and roll accordingly.

TREAT FOOT AILMENTS

This device can also be used to address Plantar Fasciitis (PF) and other foot problems. While seated, lay the textured ball on a non-slippery surface (e.g. carpet, exercise mat) and roll it back and forth beneath your foot. Focus and apply pressure on the part of the foot that hurts the most.

While the ideal therapy time is 4 to 6 minutes, you may adjust its length based on the severity or intensity of your pain. Therapy performed in the morning is most beneficial for individuals suffering from PF.

USE COLD TO COMBAT INFLAMMATION

People with mild to moderate swelling, bruising, sprain, and inflammation will benefit most from cold therapy. To freeze the balls' gel, store your Hot and Cold Massage Kit in the freezer for 4 to 6 hours. The balls will remain cold for 15 to 20 minutes. Refreeze as needed or for succeeding use.

USE HOT TO REDUCE PAIN AND RELAX MUSCLES

Heat therapy works to reduce joint pain and muscle spasms, and relax sore muscles, ligaments, and tendons. Place the balls in a container of very hot (120 degrees F or C?) tap water for about 7 to 10 minutes. Before applying directly to skin, test first to determine heat level and comfort. The balls will remain hot for about 10 to 15 minutes. Reheat as needed.

DO NOT HEAT USING OTHER IMPLEMENTS/DEVICES. The use of appliances such as conventional or microwave ovens and stove tops can result in overheating, which may in turn inflict burns on your skin and damage on your Hot and Cold Massage Kit. To avoid skin burns and damaging the seals on the balls, do not place your Hot and Cold Massage Kit in boiling water.

WARNING: Use only as directed. Discuss with your physician or physical therapist any self-treatment program you plan or wish to implement. This product is intended for people aged 13 and older. Keep out of reach of children. People with diabetes and those with open sores and cracked or blistered skin in the body part to be treated should not use the Hot and Cold Massage Kit. If symptoms get worse while using the product, reduce pressure. However, if symptoms persist, discontinue use and consult your health professional immediately.

Setup, Usage, Care, and Safety Guide

Before setting up and using your Hot and Cold Massage Kit, read carefully the following instructions.

SETUP

Thoroughly check your Hot and Cold Massage Kit to make sure that it is in excellent condition. Inspect the balls for scratches, tears, gashes, and other product defects or shipping damage. Once you confirm that both balls are in perfect shape, you may start using your Hot and Cold Massage Kit.

MASSAGE POINTERS

1. Perform each massage in a slow, even, careful, and effortless manner.
2. Do each massage as best as you can; however, do not overexert or strain yourself.
3. Breathe freely and evenly when massaging or stretching. Do not hold your breath.
4. It is best to use your Hot and Cold Massage Kit in a wide, open space where there are no furniture, appliances, and other household items that could otherwise hamper your massage.

CARE INSTRUCTIONS

1. After every use, dry your Hot and Cold Massage Kit with a clean towel. Thoroughly inspect the balls for any nicks or tears.
2. Spot clean only. Use a soft, damp cloth to remove spots, stains, or marks from your Hot and Cold Massage Kit. Air dry.
3. Do not expose your Hot and Cold Massage Kit to harsh or prickly surfaces, extreme heat, or too much sunlight. All these can damage the balls or reduce their efficacy.

SAFETY GUIDELINES


1. Use your Hot and Cold Massage Kit only as prescribed.
2. When it comes to exercise, it can never be "one type fits all." To minimize the possibility of injury, consult your physician before starting any type of exercise program or using any exercise equipment or device. This instructional sheet should in no way take the place of sound medical advice.
3. Got a heart condition, high blood pressure, or any other disease or affliction? Talk to a health professional before starting any exercise regimen or using the Hot and Cold Massage Kit.
4. Should you experience dizziness, shortness of breath, or any other discomfort while using the product, pause and rest for a moment.
5. This Hot and Cold Massage Kit does not come with any warranties or guarantees. The manufacturer and distributor cannot be held liable in the event of injury, loss, or damage resulting from the use of this product.

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