# HOT & COLD FOOT RELIEF ROLLER

Your **Hot & Cold Foot Relief Roller** is designed to soothe pain associated with Plantar Fasciitis (PF) as well as heel spurs and other common aches and discomforts. It has an ergonomic shape that conforms to the curve of your foot plus a textured surface that comfortably massages. Your roller massager is powered by a hot/cold technology, so you can expect it to stay hot or cold for as long as 20 minutes of therapy.

## Setup, Usage, Care, Massage, and Safety Guide

Prior to setup and use, please take the time to read the following guidelines:

### HOW TO SET UP

After taking it out of the box, carefully inspect your Hot & Cold Foot Relief Roller for any manufacturing defect or shipping damage such as scratches, punctures, or gashes. Once you confirm that it is in excellent condition, you may begin using your Hot & Cold Foot Relief Roller.

### HOW TO USE

#### FOR A RELAXING FOOT THERAPY

Use your Hot & Cold Foot Relief Roller to relieve PF and other foot aches and discomforts. Lay your Hot & Cold Foot Relief Roller on carpeted floor, an exercise mat, or any even, non-slippery surface. While seated, roll your massager back and forth beneath your foot. Concentrate and apply maximum pressure on those parts of your foot that ache the most. Ideally, this therapy should take anywhere from 4 to 6 minutes. However, therapy length can differ depending on pain severity. Morning treatment is highly recommended for people suffering from PF.

#### FOR A CALMING COLD THERAPY

To freeze the gel inside it, store your Hot & Cold Foot Relief Roller in the freezer for a minimum of 4 to 6 hours. While in use, it should remain cold for 15 to 20 minutes. Refreeze as needed or for future use.

#### FOR A SOOTHING HOT THERAPY

Submerge your Hot & Cold Foot Relief Roller in very hot water (about 120 degrees Celsius) for 7 to 10 minutes. To avoid burns, make sure to test the heated roller before applying directly to the area that needs treatment. While in use, it should remain hot for as long as 10 to 15 minutes. Reheat as needed or for succeeding use.

#### FOR YOUR UTMOST SAFETY AND COMFORT,

- NEVER USE OTHER APPLIANCES/DEVICES to heat your Hot & Cold Foot Relief Roller.

Using a stove top or conventional/microwave oven to heat your Hot & Cold Foot Relief Roller is extremely dangerous. It can cause overheating, which can in turn result in skin burns and damage to your roller massager. Neither should you submerge it in boiling water; doing so would destroy its seals and burn your skin.

- NEVER STAND on your Hot & Cold Foot Relief Roller. The product should only be used in a sitting position.

#### FOR A TRULY RELAXING AND EFFECTIVE MASSAGE,

- Perform every massage carefully, slowly, evenly, and effortlessly. Do not overexert.
- Do your best, but make sure you do not go beyond what your body can endure.
- Breathe freely and steadily. Holding your breath while massaging or stretching may prove harmful.
- Use your Hot & Cold Foot Relief Roller in a wide, open space. A narrow, cluttered area littered with furniture and other household items could keep you from performing and enjoying an effective massage.

#### TO KEEP YOUR ROLLER IN EXCELLENT CONDITION,

- With a clean towel, thoroughly dry your Hot & Cold Foot Relief Roller after every use.
- Spot cleaning would suffice. Using a soft, damp cloth, gently remove spots, stains, or marks from your Hot & Cold Foot Relief Roller. Air dry.
- Never expose your Hot & Cold Foot Relief Roller to harsh or prickly surfaces, extreme heat, or too much sunlight. Doing so can damage or reduce the effectiveness of your roller.

#### MORE REMINDERS TO PROTECT YOU AND YOUR ROLLER

- A particular exercise routine does not necessarily work for everybody. To avoid all forms of injury, talk to your physician before you begin any exercise program or use any exercise equipment or device. Do not substitute this instructional sheet for proper medical/professional advice.
- Do you have a heart ailment, high blood pressure, or any other disease or health issue? Make sure to see a health professional before you start an exercise regimen or use your Hot & Cold Foot Relief Roller.
- Every massage or exercise routine must be done slowly and steadily.
- If, while using your Hot & Cold Foot Relief Roller, you experience dizziness, shortness of breath, or any other discomfort, pause and rest for a few minutes.
- No warranty or guaranty accompanies your Hot & Cold Foot Relief Roller. Hence, its manufacturer and distributor cannot be held liable for any injury, loss, or damage that may result from the use of this product.

#### IMPORTANT WARNING

Use only as prescribed and described in this instructional sheet. It is best to discuss with your physician or physical therapist any at-home remedy or self-treatment method you wish to try. This product is designed for people ages 13 and older. Keep out of reach of very young children. Those with diabetes, open sores, or cracked or blistered skin in the affected body part are strongly advised against using this roller. Reduce pressure should symptoms or discomfort get worse while using the product. However, if acute symptoms continue, suspend use and call your physician immediately.

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