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How to use

The acupressure reflexology mat takes a while before getting accustomed to because of the sharp points. It may cause some discomfort or pain for few minutes until your feet becomes stable. With its pointed nibs, it massages your unique nerve endings.

1. Step on the acupressure reflexology mat for 30 to 60 seconds or until you are able to stay longer.
2. Change your position from side to side using your weight to emphasize the foot area you want to stimulate.
3. Repositioning your feet also helps initiate movement and alleviates pain associated with plantar fasciitis.
4. If the pain becomes too much to handle, please step back until you feel comfortable enough to step in once again.
5. You can also try sitting while pressing your feet on the acupressure reflexology mat. That way you can manage the pain easily if its too much to handle standing with full weight pressing down.

The acupressure reflexology mat relaxes the feet muscles helps in reducing stress without causing any side effects, designed for all foot types it helps relieve conditions such as plantar fasciitis, heel spurs, arch pain, myofascial pain syndrome, neuropathy including tired and sore feet from standing too long or due to activities.

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May not be suitable for people with underlying health conditions. Please discontinue use if you feel extreme pain. Consult a physician prior to engaging in physical activity.

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